

# RAPYDTEST®

FOR THE DETECTION OF H.PYLORI AG IN FAECES



H.pylori Ag  
RAPYDTEST®

## Intended Use

The H.pylori Ag Rapydtest® is a lateral flow chromatographic immunoassay for the qualitative detection of H.pylori antigen in human faecal specimen. It is intended to be used by professionals as a screening test and as an aid in the diagnosis of infection with H.pylori.




## Performance Characteristics

### Clinical Performance

324 faecal samples collected from subjects with symptomatic gastrointestinal disorders and non-gastrointestinal symptoms were tested with the H. pylori Ag Rapydtest® and with the UBT as reference test.

A comparison of the results for all subjects is shown in the table below.



UBT	H.pylori Ag Rapydtest®		TOTAL
	POSITIVE	NEGATIVE	
POSITIVE	118	7	125
NEGATIVE	0	199	199
TOTAL	118	206	324

Relative Sensitivity: 94.4%

Relative Specificity: 100.0%

Overall Agreement: 97.8%



**BACTERIOLOGY**

SINGLE USE IN VITRO DIAGNOSTIC DEVICE



H.pylori Ag  
RAPYDTEST®


### Reagents and Materials Provided

1. Individually sealed foil pouches containing:
  - a. One cassette test device.
  - b. One desiccant.
2. Sample extraction tubes, each containing 2ml of extraction buffer.
3. Plastic droppers for transferring watery stool.
4. One package insert (instruction for use).


### Specimen Collection and Handling

Consider any materials of human origin as infectious and handle them using standard biosafety procedures.

### Procedure for Stool Sample Collection



### Test Procedure



**Do not read results after 20 minutes. To avoid confusion, discard the test device after interpreting the result.**

### Interpretation of Assay Result

1. Negative Result:  
If only the C band is developed, the test indicates that no detectable H.pylori antigen is present in the specimen.  
The result is negative.



2. Positive Result:  
If both C and T bands are developed, the test indicates the presence of H.pylori antigen in the specimen.  
The result is positive.



3. Invalid:  
If no C band is developed, the assay is invalid regardless of any colour development on the T band as indicated below. Repeat the assay with a new test device. Excess faecal specimen can lead to invalid test results; if this is the cause, re-sample and re-test (see instructions for collection of specimen).



### References

1. Marshall BJ, et al, pyloric campylobacter infection and gastroduodenal disease. Med.J. Aust. 1985, 142: 439-444
2. Lambert IR, Lin SK, and Aranda-Michel J, helicobacter pylori Scan. J. Gastroenterol. 1995, 30 suppl 208: 33-46
3. Vans DJ, Evans DG, et al A sensitive and specific serologic test for detection of campylobacter pylori infection. Gastroenterology. 1989, 96:1004
4. Shimoyama T, Kato C, et al. Applicability of a monoclonal antibody- based stool antigen test to evaluate the results of Helicobacter pylori eradication therapy. 2009, May 62(3): 225
5. Krause R, Muller G, Doniec M. Evaluation of a rapid new stool antigen test for diagnosis of Helicobacter pylori infection in adult 1008, 46(6): 2062
6. Altman E, Fernandez H. et al Analysis of Helicobacter pylori isolates from Chile: Occurrence of selective type I Lewis b antigen expression in lipopolysaccharide. 2008, 57(pt 5): 585

PRODUCT	PACK SIZE	CODE
H.pylori Ag Rapydtest®	25	1632

Products can be ordered direct from Apacor or from an appointed distributor  
Visit our website for all the latest information [www.apacor.com](http://www.apacor.com) or email on: [sales@apacor.com](mailto:sales@apacor.com)



UNIT 5 SAPPHIRE CENTRE  
FISHPONDS ROAD, WOKINGHAM  
BERKSHIRE, RG41 2QL, ENGLAND  
TEL: +44 (0)118 979 5566  
FAX: +44 (0)118 979 5186

