




PARASITOLOGY



Protozoa and cysts found in faeces




Helminth eggs and larvae found in faeces




Stains


Stains for faecal deposits: Lugols Iodine



Staining of faecal smears: Rapid Field stain, Iron haematoxylin or Giemsa stains




Modified Trichrome stain for Microsporidia




E intestinalis - Showing vacuole and suture


Modified Zn for Cryptosporidia



Cryptosporidia



Midi Parasep®
Faecal Concentrator



Mini Parasep® SF
Solvent Free System

Products can be ordered direct from Apacor or from an appointed distributor
Visit our website for all the latest information www.apacor.com or e-mail on: sales@apacor.com